

Diagnostics and Imaging Equipment



Treatment Units



Handpieces and Instruments



Endodontic Systems



Laser Equipment



Laboratory Devices



Signo G10 II

Creating new freedom in dentistry



Distributed by
J. Morita Europe GmbH
Justus-von-Liebig-Str. 27a
63128 Dietzenbach
Germany
T +49. 6074. 836 0
F +49. 6074. 836 299
www.morita.com/europe

Developed and manufactured by
J. Morita Tokyo Mfg. Corp.
7129 Komuro Ina-machi,
Kitaadachi-gun, Saitama
362-0806 Japan

Subject to technical changes and errors.
JME EN FA 0815 *2.5

Signo G10 II – different versions for individual requirements



Double growth for the Morita Signo series which has proven its worth time again and is extremely popular among both dentists and dental assistants alike: Convincing features of the new Signo G10 II treatment unit include not only its appealing design but also, and above all, its excellent ease of maintenance, multiple functions and tremendous versatility. Both designs score with a wide range of high-quality instruments – from turbines to scalers up to micro motors, all of which are fitted with LED lights. In addition, the high-powered, long-life and efficient “LunaVue EL II” LED OP light is a guarantee for a perfect view of the treatment area at all times.

A particular highlight of the whip-arm delivery is the operator’s element with 120 degree swing movement without any irritating ratchet function. Instrument balance is optimized by the DUAL BALANCE system with its dynamic whip-arm strain relief. In addition, the variable deposit tray means that the operator’s element can be placed in a variety of positions so that it can cover all usual treatment positions for the practitioner. The over-the-patient delivery also offers a large range of benefits – including the over-the-patient arm with a swing range that has been improved by 20 degrees.




Whip-arm delivery

- 120° swing range for very large freedom of movement for treatment
- Long instrument hoses for high activity radius
- Flexible setting with various deposit trays as storage area for instruments
- Side handles or optional central handle for easy positioning



Over-the-patient delivery

- Can be fitted with up to 5 instruments
- Swing range optimized by 20° for flexible positioning
- Easy access to instruments



Easy patient positioning

- Hydraulic drive for gentle patient chair movement
- Exact focus on the treatment area
- Automatic height adjustment of cuspidor and OP light



Intuitive operation

- Pushbutton control of the patient chair and instrument parameters
- Easily comprehensible symbols for improved overview
- As an alternative, the main parameters can be controlled by foot control
- Up to 5 different treatment programs can be stored individually



Comfortable upholstery

- Soft upholstery for advanced patient comfort
- Headrest adjustable to the size of the patient for comfortable and stable position of the patient’s head
- Optional standard or multi-jointed headrest